

	TURNOUT 2	TURNOUT 4	OPTIONS
Height x width x length:	38x450x480 mm	60x450x510 mm	<ul style="list-style-type: none"> <li>● BEV-seat</li> <li>● BEV-seat accessories</li> <li>● Rails for BEV-seat.</li> <li>● Vehicle model specific mounting brackets</li> <li>● Universal mounting bracket</li> <li>● Remote control (wireless)</li> <li>● Tilda, the seat tilter</li> </ul>
Seat level above threshold, incl. BEV-seat and rails:	90 mm	90 mm	
Maximum load Turnout	250 kg	250 kg	
Maximum load Turnout E	150 kg	150 kg	
Weight, model 2 and 4	18 kg	20 kg	
Weight, model 2E and 4E	18,8 kg	20,7 kg	

Fulfills following regulations: I74/408/EEG to 96/37/EG, 93/42/EEG, ISO 3795 FMVSS/CMVSS 207, 208, 210 and 302 (Turnout with BEV-seat)

### THE SWIVEL BASE

Turnout is designed to be easy to use, strong and dependable. It is made out of 8 mm coated Domex 355 steel and has a plastic cover at the rear edge. The swivel has a swedish made, large dimensional axial ball bearing to endure strong force.




### THE SEAT

Turnout is combined with a low profile seat, such as the BEV seat. The BEV seat has many accessories that improve the seating comfort and support for the user.

### MOUNTING BRACKET

A mounting bracket, adapted for the specific vehicle model, is fixed in the original seat mountings in the floor.



# TURNOUT®

Swivels out the seat and makes it easier entering and exiting the vehicle



Please find more information on [www.autoadapt.com](http://www.autoadapt.com)



**Autoadapt AB**, Åkerivägen 7,  
S-443 61 Stenkullen, Sweden  
Tel: +46 302 25400 Fax: +46 302 25407  
E mail: [contact@autoadapt.se](mailto:contact@autoadapt.se) [www.autoadapt.com](http://www.autoadapt.com)

**Autoadapt UK Ltd**, Unit 1, Windsor Industrial Estate,  
Rupert Street, Aston, B7 4PR Birmingham, UK  
Phone: 0121 333 5170 Fax: 0121 359 6194  
E mail: [contact@autoadapt.co.uk](mailto:contact@autoadapt.co.uk) [www.autoadapt.co.uk](http://www.autoadapt.co.uk)



INVENTIONS FOR QUALITY OF LIFE

# TURNOUT® - simplifies getting in and out of the vehicle

## Manually operated or powered

The Turnout swivel base is of great help for those who have difficulties getting in and out of their vehicle. By swivelling the seat outwards you have considerably improved your possibilities of sitting down or rising up.

**TURNOUT 2** (2-door vehicles) and **4** (4-door vehicles) are manually operated. Lift the lever, placed at the side of the seat, to release it and then swivel the seat outwards.

**TURNOUT 2E** and **4E** are powered and operates by pressing a hand control button. The big advantage of having a powered Turnout is that you can fully concentrate on moving your feet over the sill. Just let go of the button if you need to interrupt the swivelling movement.



## Low profile vehicle seat

To preserve a suitable seat height inside the vehicle the Turnout shall be combined with a low profile seat, such as the BEV-seat. It is not recommendable to use the original seat with the Turnout.

## No damages to the vehicle

The Turnout is mounted in the original seat bolt holes in the floor. When the vehicle needs to be replaced the original seat can be mounted again and the vehicle has not lost resale value.



## Turnout combined with Tilda, tilting the seat in three steps

A TURNOUT combined with TILDA makes it even easier to get in and out of the vehicle. Tilda provides the possibility to tilt the entire seat backwards in three steps. By resting the feet on the optional footrest, they will comfortably pass over the vehicle threshold while swivelling inside. Tilda is crash tested for sitting in tilt position one during drive.

## Safety

The Turnout swivel base is crash tested and meets current safety regulations. The product is CE-marked.



## Combined with Carony

Turnout is also an integral part of the Carony system, which is a wheelchair with a seat that can be moved on rails from the wheel unit to the vehicle swivel base.

By using a Carony the wheelchair occupant does not need to be lifted in and out of the vehicle. The user can be seated in the same seat from inside the house and right into the vehicle front seat position.

The Carony wheelchair is easy to manoeuvre and functions as an ordinary comfort wheelchair. It is available with 24" driving wheels or as a transportation wheelchair with 12" rear wheels.



*The Carony Kids is a Carony version which adapts to the child as it grows.*



*Carony 24" has driving wheels that can easily be removed. The unit and the wheels can then be stored in the vehicle luggage area.*